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



Tootja

REF Referents = toote number

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Operation Instructions: Blood gas capillaries

100µL, 105µL, 130µL, 100µL, 105µL, 130µL 140µL, 165µL,
185µL, 200µL, 235µL, 280µL (Sodium heparinised)

Instructions for use

The capillaries are disposable glass capillary pipettes intended for one-time use with subsequent disposal.

Suitable for medical professionals only.

Application temperature: +2 – +40°C

National guidelines for occupational safety and quality assurance must be observed.

Purpose

✓ **Intended purpose:** Blood gas capillaries Na-Hep for transferring defined volumes of capillary whole blood for quantitative analysis. The heparin contained in the glass capillary prevents the coagulation of blood in order to avoid a negative impact on the downstream analysis. The blood is entered into the capillary directly at the analysis device.

✗ **Prohibitions of use:** Blood gas capillaries, Na-Hep are NOT intended for the following applications:

Sampling of other liquids than capillary whole blood, storage of samples, measurements within the capillary itself.

Hazards and safety

The precautionary measures required for the use of laboratory reagents and body fluids as well as microbiological samples must be observed. Handling may only be carried out by qualified personnel. National and internal laboratory guidelines for occupational safety and infection control must be observed. Suitable protective clothing, protective goggles and disposable gloves must be worn during work.

Attention must be paid to effective protection against infection in accordance with the laboratory guidelines.

The user is committed to report any serious incidents in connection with the product to the competent authority of the Member State in which the user and/or the patient is registered and to Hirschmann Laborgeräte GmbH & Co. KG.

The capillaries comply with the requirements of Regulation (EU) 2017/746 regarding in vitro diagnostic medical devices.

Disposal

The national statutory regulations must always be observed when disposing of the capillaries. Used capillaries contaminated with potentially infectious substances of human origin may pose a risk of injuries, infections or microbiological hazards. Potentially contaminated capillaries must be disposed of in such a way that injuries, infections or microbiological hazards to third parties are excluded. Unused capillaries can be disposed of with the household waste. Unused capillaries must be disposed of in such a way that injuries to third parties are excluded.

No.	Description
9221210	BLOOD GAS CAP. NA-HEP 130 µl; L 100MM; OD 1,75MM; ID 1,35MM
9221212	BLOOD GAS CAP. NA-HEP 165 µl; L 125MM; OD 1,75MM; ID 1,35MM
9221215	BLOOD GAS CAP. NA-HEP 200 µl; L 150MM; OD 1,75MM; ID 1,35MM
9221275	BLOOD GAS CAP. NA-HEP 100 µl; L 75MM; OD 1,75MM; ID 1,35MM
9221280	BLOOD GAS CAP. NA-HEP 105 µl; L 80MM; OD 1,75MM; ID 1,35MM
9222210	BLOOD GAS CAP. NA-HEP 185 µl; L 100MM; OD 2,35MM; ID 1,55MM
9222212	BLOOD GAS CAP. NA-HEP 235 µl; L 125MM; OD 2,35MM; ID 1,55MM
9222215	BLOOD GAS CAP. NA-HEP 280 µl; L 150MM; OD 2,35MM; ID 1,55MM
9222275	BLOOD GAS CAP. NA-HEP 140 µl; L 75MM; OD 2,35MM; ID 1,55MM

IVD In vitro diagnostic medical devices

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