HIRSCHMANN®

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Instruction manual

Please read before use!

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Contents

	Page
1. Warning and safety instructions.....	5
2. Individual components	9
3. Function overview.....	10
4. Prohibited use	12
5. Restrictions on use.....	13
6. Preparing the unit.....	14
6.1 Using the UV protection window.....	15
6.2 Securing the suction hose	16
6.3 Securing the recirculation hose	17
6.4 Securing the ejection unit	17
6.5 Screwing the unit onto the bottle.....	18
6.6 Turning the ejection unit to the correct position.....	19
6.7 Bottle label aligning unit	19
6.8 Battery installation (optional).....	20
6.9 Battery charging (optional)	21
7. Start of work	22
7.1 Overview of symbols	23
7.2 Activating/deactivating the unit.....	24
7.3 LED status	26
7.4 Modification of decimal places (10 ml)	26
7.4.1 Modify displayed decimal places.....	27
7.5 Closure cap removal.....	27
7.6 Venting the unit	28
8. Titration	29
8.1 Filling	30
8.2 Configuring of suctioning and ejection speed	32

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	Page
8.3	Rapid titration..... 33
8.4	Droplet titration..... 35
8.5	Clear..... 36
9.	Special functions 37
9.1	Main screen..... 37
9.2	Language selection..... 38
9.3	Quick-Cal..... 38
9.4	PIN Protection 40
9.4.1	Activate PIN protection..... 40
9.4.2	Deactivation of the PIN..... 41
9.4.3	Activation of restricted user PIN 41
9.4.4	Deactivation of restricted user mode..... 43
9.4.5	If the PIN is forgotten 43
9.5	Method saving..... 44
9.6	Activating the method..... 45
9.7	Switches 46
9.8	Info..... 47
11.	Interface 48
11.	Cleaning..... 48
11.1	Service/Cleaning 50
11.2	Basic cleaning with the venting function 51
11.3	Intensive cleaning..... 52
12.	Sterilisation..... 55
12.1	Preparation..... 56
12.2	Dismantling 56
12.3	Loosening the suction, ejection and recirculation valve 57

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	Page
12.4 Sterilisation procedure	60
12.5 Assembly	61
12.6 Tightening the suction, ejection and recirculation valve.....	62
12.7 Function inspection.....	63
13. Valve replacement	64
13.1 Suction valve replacement	64
13.2 Ejection valve replacement	65
13.3 Recirculation valve replacement	65
14. Repair by manufacturer	66
15. Volumetric inspection.....	66
16. Trouble-shooting	72
16.1 Messages.....	72
16.2 Other problems - solutions.....	77
16.3 Contacting the manufacturer.....	80
17. Technical data.....	83
18. Material	85
19. Accessories and spare parts	86
20. Clearance certificate form.....	88

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1. Warning and safety instructions

Warning!

Read the instruction manual prior to using the unit and observe all safety instructions. Any person working in the laboratory should be familiar with these safety instructions prior to using the unit, and the safety instructions should be accessible at all times. However, the instruction manual cannot illustrate every safety problem which may occur. Every user is personally responsible for observing the safety and health protection regulations and determining restrictions prior to using the unit. These restrictions should then be observed.

- Observe all general safety instructions for the laboratory e. g. wearing of protective clothing, safety glasses and protective gloves when using certain liquids.

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1. Warning and safety instructions

- Observe the instructions and regulations of reagent manufacturers.
- The unit should only be used for titrating liquids, and it is imperative that usage prohibitions and restrictions be observed (consult the manufacturer if necessary).
- Prior to use, check all connections for leaks and ensure that the ejection unit and suction hose are securely fitted. A loosely-fitted ejection unit could cause splashing when the piston unit is lowered.
- It is imperative that the closure cap should be removed from the ejection hose prior to titrating. Never activate the unit while the closure cap is fitted.
- Ensure that neither you yourself nor other persons are at risk. Never work facing in the direction of another person. Avoid splashing and use suitable collecting vessels. Always titrate onto the interior wall of the collecting vessel rather than directly onto its base or the contents contained therein.

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1. Warning and safety instructions

- Never use force!
- Seal the ejection hose with the closure cap after use. **Caution! Reagent can drop out!**
- Use only original accessories and original spare parts.
- Allow components to cool to room temperature after sterilisation (see point 12).
- In case of malfunction, e. g.
 - piston unit stiffness or seizure
 - sticky or leaking suction, ejection or recirculation valve
 - leaks
 - component breakagestop working immediately. Observe the instructions on cleaning (see point 11) and problem solutions (see point 16) before continuing to use the unit. Send to manufacturer for repair if necessary. Observe point 14 in this respect!

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
1. Warning and safety instructions

- Do not deposit the unit or power supply unit in water or other liquids to avoid the risk of electric shock!
- Ensure that the mains power voltage corresponds to the details on the power supply unit.
- The power supply unit should not be used any further if the cable or power supply unit are damaged in any way.
- Never touch the plug on the power supply unit with wet hands!
- The power supply unit should only be used indoor. Protect against damp!

Warning!

Liability is excluded for the unit or any subsequent damage if the unit is opened. Avoid incorrect use! Always use an original power supply unit and original accessories and spare parts to avoid damage to the unit.

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

2. Individual components

1. Ejection unit (complete)
2. Recirculation valve
3. Suction valve
4. Ejection valve
5. Basic unit
6. Control unit
7. Power supply unit

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3. Function overview

- The unit is a motorised liquid titrating system. The basic unit (A 45 thread) and adapter variants (see accessories) facilitate the use of commercially available reagent bottles.
- The system employs the piston stroke principle and functions via ball valves. The piston unit is advanced by a DC motor with gear drive and increments. Actuation is realized with a microprocessor control. Power is supplied by a power supply unit. The unit can be optionally powered with a battery pack.

Note:

The upper part of the basic unit is not sealed against gas and damp. It is therefore not permissible to operate the unit in caustic fumes.

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3. Function overview

- The suction valve is opened when the piston unit is drawn upwards. Liquid is suctioned into the cylinder from a storage bottle through the suction hose. The suction valve closes when suctioning is completed.
- The ejection valve is opened when the piston unit is lowered and the liquid in the cylinder is ejected.

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4. Prohibited use

The unit should not be used with the following:

- Liquids which attack the materials Halar (ECTFE), FEP, PTFE, borosilicate glass 3.3 or Hastelloy.
 - Suspensions (e.g. active carbon suspension), as solid particles can block the valves.
 - Hydrofluoric acid, as it attacks glass.
 - Combustible media (note the flashpoint).
 - Use in potentially-explosive environments.
-

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5. Restrictions on use

- The working temperature range for the unit and liquid lies between 4 °C and 40 °C.
- The unit should only be used and optional batteries loaded in areas where it can be assured that no hazardous potentially explosive atmospheres occur.
- Do not deposit in potentially explosive areas.

Attention - important information!

Never use force if the piston unit has seized. Terminate the process immediately and realize basic cleaning (see point 11).

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6. Preparing the unit

Warning, please note before assembly!

- Check all parts for evidence of damage!
 - Flawless and safe operation of the unit is only assured if the undamaged original accessories included in the delivery are used.
 - Never use damaged accessories or accessories gloves.
-

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6.1 Using the light protection window

The existing clear inspection window can be replaced with the coloured UV protection window to protect media which are sensitive to light.


The clear inspection window should be raised initially on one side for removal to ensure that it slides out of its fixture, and then raised again on the other side.

The coloured UV protection window can now be suitably fitted (note correct position).

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
6.2 Securing the suction tube


- Hold the suction tube beside the storage bottle. The hose end with the union nut should face towards the base of the bottle.
- Cut the hose at a slant level to the lower edge of the locking ring on the storage bottle.

- Push the suction hose on up to the stop.
- Then tighten the union nut onto the suction valve.

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

6.3 Securing the recirculation hose

Insert the recirculation hose firmly into the recirculation opening in the valve block.

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
6.4 Securing the ejection unit

- Ejection unit with:
 - hose guide
 - union nut
 - ejection hose
 - closure cap
- Push the ejection unit on up the stop.

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6.4 Securing the ejection unit

- Tighten the union nut firmly and ensure that it is securely fitted.

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
Attention!

- Ensure that the ejection unit is securely fitted.
- Only use original ejection unit (see point 19, Accessories and spare parts).
- Do not connect a deformed or damaged ejection unit.

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6.5 Screwing the unit onto the bottle

Screw the unit onto the storage bottle (use the threaded adapter where appropriate!)

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6.6 Turning the ejection unit to the correct position

The unit and ejection unit can be aligned in practically any direction relative to each other, thanks to two rotation axes integrated in the unit.

Warning!

- Observe all safety instructions, usage prohibitions and restrictions on use.
- Wear protective clothing, safety glasses and protective gloves.

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6.7 Bottle label aligning unit

- Turning the unit into the desired position relative to the bottle label.
- The unit should not be positioned directly at the right-hand stop of the rotation axis (see point 12.2 for dismantling position).

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6.8 Battery installation (optional)

The control unit should be opened to facilitate installation of the battery pack. All cables should be removed from the control unit prior to commencing battery installation. All 4 screws at the rear of the control unit should then be loosened (Phillips screw PZ1).

The control unit should then be rested on its feet again. The electronic should be tilted away over the upper edge of the housing after removing the cover. Attention should be paid to the wiring during this. The battery pack is then fixed in the centre of the housing floor with the assistance of the double-sided adhesive tapes provided.

Ensure that the connection cables for the battery are in the upper right-hand corner during this.

The battery connection plug should be connected to the respective connection on the PCB.

The PCB can then be positioned on its base, the cover replaced and screws hand-tightened again. Attention should once more be paid to the wiring during this.

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6.9 Battery charging (optional)

The unit is, in principle, ready for immediate use when connected to the mains power supply. However, the battery should be charged prior to initial commissioning if operation with the battery option is chosen.


The battery is charged when the power supply unit is connected, regardless of whether the unit is activated or deactivated.

Charging is indicated by the battery indicator on the display and by the flashing red LED on the display module if the unit is switched off. Charging is terminated automatically when the battery is fully charged (battery overloading is impossible!) and a green LED flashes on the display module. The unit is also ready for operation during battery charging.

A warning message appears on the display when the battery capacity is no longer adequate to ensure correct operation. Work may only be continued in this case after the power supply unit is connected.

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
Warning:

The unit should only be operated with the original power supply unit. The mains voltage must correspond to the voltage specified on the power supply unit! (Max. deviation $\pm 10\%$).

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7. Start of work

Connect the basic unit and control unit with the cable provided (blue 5-pole connector and socket).

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Establish the power supply by connecting the power supply unit (or, optionally, the battery pack). The unit is now ready for operation.

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Not all buttons can be used when certain operating conditions prevail, particularly if a process (titration, venting, etc.) is currently active. These are then temporarily deactivated. Inactive buttons are illustrated with a grey symbol margin.

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7.2 Activating/deactivating the unit

The unit is activated by pressing the touch screen for approx. 2 s. The following screen then appears:

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To switch off the unit, activate the menu , and select the “Off” menu item.

Acknowledge with .

- If the unit is only supplied with power via the optional battery, it deactivates automatically after 5 minutes in the absence of entries by the user or if no process is currently active.
- All user configurations are retained if the unit is deactivated automatically or with the “Off” button, so that work can be continued immediately by pressing the touch screen again.

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Warning!

- Safety instructions must be observed (see point 1)!
 - Usage prohibitions and restrictions of use must be observed (see point 4 and 5).
 - The unit should only be utilised with the ejection unit turned away from the user and other persons. The ejection unit and valve block should be positioned turned away to the left or right of the body for this purpose.
-

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7.3 LED status

Basic unit	
LED flashes green	Unit ready for operation
LED flashes red	Error (see point 16)

Control unit	
LED flashes green	Power supply unit plugged in or power supply unit plugged in and battery fully charged (only when the unit is deactivated)
LED flashes red	Battery charging (only when the unit is deactivated)

7.4 Modification of decimal places (10 ml)

In addition to the 50 ml and 20 ml unit the 10 ml unit allows slower velocities and greater accuracy to be achieved.

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7.4.1 Modify displayed decimal places

The volume for the 10 ml unit can be displayed to either 2 or 3 decimal places. To switch from one to the other call up the menu  and press the button .

The main display with the changed setting will then appear. When switching from 3 decimal places to 2 places the third place will be cut off.

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7.5 Closure cap removal


- Position a collecting vessel under the ejection unit.
- Pull of the closure cap as shown in the illustration. **Warning: reagent can drop out during this operation!**

- Push the closure cap to the rear as shown in the illustration. The closure cap can remain in the retention position while work is being realized.

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7.6 Venting the unit

The recirculation lever should be positioned at
recirculation (the tip of the recirculation lever
points backwards) if air bubbles visibly form in
the cylinder chamber (e. g. after screwing onto
a new storage vessel and subsequent filling).

The unit is then automatically vented using the “Venting”
symbol in the main window.

1 ml is suctioned and the cylinder drained when the
symbol is pressed. 40% of the nominal volume is then
suctioned and dispensed again. Then “Venting” symbol
is converted into a “STOP”  symbol during venting.
The venting process can be interrupted by pressing the
“STOP” symbol.

Venting is not compulsory if the unit was correctly deacti-
vated after its last use and the control unit not detached
from the basic unit. However, venting is recommended in
this case.

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7.6 Venting the unit

It is imperative in some cases that venting will be realized first when activating the system. This is indicated in that only the “Venting” and “Menu”  buttons are then active.

8. Titration

During titration, the volume dispensed through the different functions is continually added. This value can be viewed on the display when the titration end point is reached.

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8.1 Filling

- Pressing the “Filling” symbol  in the main window fills the glass cylinder with the nominal volume of 10, 20, or 50 ml.
- It is possible to fill the cylinder with a randomly-selected volume. The volume to be suctioned is configured in the main window prior to filling for the purpose.

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All menus are laid out with the same schematic:

- An input field opens when the volume indication is pressed. The current value is shown immediately beneath. The minimum and maximum value are displayed adjacent to this, and these values should not be exceeded or undershot.
- Inputting of the new value is realized via the numerical sequences, and incorrect entries can be corrected with the “DEL” field.

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8.1 Filling

- The new value is imported by pressing the “OK” ✓ symbol, provided it is valid. Please note the min. and max. values in the case of an invalid value.
- The new value is rejected by pressing the “Abort” ↶ symbol and the main window opens.

If the entry is acknowledged with “OK” ✓, the volume now configured is displayed with an “*” as indicator to the right of the display. This indicator illustrates that not the titration volume, but rather the preselected volume is currently shown on the display. The cylinder is filled with the preselected volume using the “Filling” ⬆️ symbol in the main window.

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8.2 Configuring the suctioning and ejection speed

It is possible to configure the suctioning and ejection speed independently of each other. The current values are illustrated in (ml/s) on the right of the main screen beside the respective symbols. The speeds configured can vary slightly depending on the ambient and operating conditions.

- An input field is opened in which the parameter is shown at the top when the symbol for the suctioning speed  and ejection speed  is pressed. The current value is shown immediately beneath. Details of the minimum and maximum value are displayed adjacent to this, and these values should not be exceeded or undershot.
- Inputting of the new value is realized via the numerical sequencs, and incorrect entries can be corrected with the “DEL” field.
- The new value is imported by pressing the “OK”  symbole, provided it is valid. Please note the min. and max. values in the case of an invalid value.
- The new value is rejected by pressing the “Abort”  symbol and the main window opens.

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8.3 Rapid titration

- Rapid titration without volume preselection: The titration volume is dispensed by pressing the  symbol. The volume value currently dispensed is indicated in the volume line. The process is stopped immediately if the button is released.
- Rapid titration with volume preselection: A volume previously entered via the numerical sequence is dispensed after the rapid titration  button is pressed. The volume is configured here as described in point 8.1.
- If the entry is acknowledged with “OK” , the volume now configured is displayed with an “*“ as indicator to the right of the display. This indicator illustrates that not the titration volume, but rather the preselected volume is currently shown on the display.

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Note

The system automatically generates the last 10% of the preselected volume in droplets. The value dispensed is displayed in the volume line on completion.

Important user advantage:

The STOP button  can be used during titration with a preselected volume to stop the current volume dispensing immediately. The volume value dispensed is shown in the volume line. Further working is, of course, possible with the “Rapid titration”  or “Droplet titration” button  at any time.

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8.4 Droplet titration

Pressing the “Droplet titration”  symbol in the main window dispenses the volume in droplets. Droplet dispensing is accelerated if the button is held for a longer period. The titration volume already dispensed can be read continually on the volume display. The smallest titration volume indicated is 0.01 ml. The process is stopped immediately if the button is released.

Note

It is sensible to change to droplet titration  before the colour change is reached to avoid any excessive titration.

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Note

If the liquid still contained in the glass cylinder is not adequate for the configured volume, the cylinder is automatically refilled. The available addition function continually adds the volume already dispensed further.

8.5 Clear

After pressing “CLEAR” , the display for the volume dispensed is reset to the start setting (0.000 ml). A new titration procedure can thus be started without previously emptying or refilling the cylinder anew.

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9. Special functions

These can be activated with the “Menu”  button.

9.1 Main screen

The unit is set to default mode by pressing the “Default settings” button. All user-specified settings are deleted. Resetting can be acknowledged or interrupted in the subsequent message. You can return to the main screen via the “Main” menu item.

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
9.2 Language selection

The user can choose from four different display languages:

- German
- English
- French
- Spanish


Configuration of the display language is realized in the “Language” submenu.

The unit imports the new language after one of the language fields is pressed.

9.3 Quick-Cal

The unit is configured by the manufacturer for a reference temperature of 20°C, taking the cubic thermal expansion coefficient for borosilicate glass (DIN ISO 4787) into consideration (EX 20°C). Quality 3 water is used for this adjustment, conforming to DIN ISO 3696.

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
Customer calibration enables the user to calibrate his unit himself to exactly suit a special application case. 10 ml, 20 ml or 50 ml (depending on the unit) of the medium is dispensed for this purpose and then weighed on an adequately-accurate scale, with actual volume being calculated taking the density into consideration. The actual volume can now be entered in the

“CAL” submenu. It is important to note here that the value calibrated by the customer should not deviate more than 10% from the nominal value (an error message is generated if the deviation is too high). After successful activation of the value calibrated by the customer by pressing the “CAL” button, the right-hand display changes to “ON”, and the “CAL” indication also appears in the main window adjacent to the volume. The unit now automatically corrects each volume to be dispensed on the basis of the calibration realized. The value calibrated by the customer can also be deactivated again pressing the “CAL” button again, and the unit now continues working with the default setting. The calibration value entered once is retained after the unit is switched off and can be activated or deactivated at any time.

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
9.4 PIN Protection

All units are provided with a function to restrict user access. If entry of a PIN is required after turning the unit on then this PIN can either prevent the entire unit from being used, or prevent settings from being changed by unauthorized persons. This function is deactivated in the factory setting.

9.4.1 Activate PIN protection

To activate the PIN protection, select the menu  and then press the button “Administrator PIN” . A four-digit PIN can now be entered and confirmed by pressing the  button.

The next time the unit is switched on the same display will appear. If the PIN is entered correctly and confirmed by pressing the tick button then the unit can be used and all functions will be available. If the correct PIN is not entered then it will not be possible to operate the unit.

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9.4.2 Deactivation of the PIN

If the unit protected through a PIN then this may be deactivated using the Administrator PIN button. The code “0000” is then entered and confirmed by pressing the “OK” button. This fully deactivates the PIN protection.

9.4.3 Activation of restricted user PIN

Once an Administrator PIN has been assigned an User PIN can be additionally assigned through the button . Anyone who enters this PIN when using the unit will only be able to access limited functions.

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9.4.3 Activation of restricted user PIN

Only the following functions can be accessed for the titration unit:

- Rapid titration
- Dropwise titration
- Clear
- Fill
- Adjust volume
- Venting
- Menu → Switch off
- Menu → Info

If the Administrator PIN and User PIN are the same then the unit will react as if no User PIN has been assigned.

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9.4.4 Deactivation of restricted user mode

To deactivate an User PIN already assigned the User PIN button is pressed. It is then necessary to enter the code „0000“ in the display and confirm this. This will fully deactivate the restriction. If the administrator PIN is deactivated as described at point 9.4.2 above then the User PIN will be automatically removed as well!

9.4.5 If the PIN is forgotten

If the User PIN has been forgotten then the user can log in using the Administrator PIN and a new PIN can be entered using the User PIN Button. If the Administrator PIN has been forgotten then it is necessary to send the serial number of the dispenser stamped on its base to Hirschmann Laborgeräte, for instance by emailing it to info@hirschmannlab.com. A new PIN will then be sent to you.

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9.5 Method saving

It is possible to save all settings relating to volumes (or partial volumes), current customer value calibration, active or inactive customer value calibration, suctioning and ejection speed as a method. Nine memory locations are available for this purpose.

Press “Save method” in the “Method” submenu to save a method. A numerical field appears in which you can select a memory location (1-9) in which the method should be saved. Memory locations which are already occupied are indicated with red numbers, but these can be written over. The memory location selected appears above the numerical field. After pressing “OK” ✓ the method is saved to the selected location. The indication M1-9 (depending on the saved method) now appears in the main window to the left of the volume indication.

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9.6 Activating the method

In order to activate a saved method, press “Activate method” in the “Method” submenu. A numerical field appears in which the method to be activated can be selected. Unoccupied memory locations have a grey background and cannot be selected. Pressing “OK” ✓ imports the new method and you change to the main window. The indication M1-9 now appears in the main window to the left of the volume indication.

Note

The M1-9 indication disappears on the display again as soon as a value has been changed after saving or activation of a method, as work is no longer realized according to this method. The same applies if work is realized alternately with a 10, 20 and 50 ml basic unit using the same control unit. This means that each unit type has its own 9 memory locations. For example, a method saved with a 20 ml unit cannot be activated with a 50 ml unit.

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9.7 Switches

The connection of additional devices enables external actuation of certain button on the main screen. The following devices are currently available for this purpose:

- Pedal switch
- Hand switch
- Mouse with two buttons
- Control cable

The auxiliary device is connected to the control unit. Only one auxiliary device can be connected to each control unit. The pedal and hand switch each have one button function, while the mouse and control cable can be configured with two button functions.

In order to configure this function, the auxiliary device should be first connected to the black interface provided on the control unit. Then press the “key 1” or “key 2” button in the “key” submenu. These buttons contain a list of file auxiliary devices that can be configured with the “key 1” or “key 2” buttons. A further submenu opens in which the function can be assigned to the button by pressing one of the symbols.

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9.7 Switches

The selection now appears with a red margin. The procedure should then be acknowledged with “OK” ✓.

The following functions can be assigned during titration:

- Rapid titration
- Droplet titration
- Clear
- Filling
- Venting

The assigned button function can be actuated from the external unit and the touch screen as soon as an external device is connected.

9.8 Info

Serial and revision numbers for both the control unit and basis unit can be viewed in the “Info” submenu.

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10. Interface

The unit can be operated from a PC via a serial or USB interface. A separate cable is available for this purpose which should be connected to the red interface at the rear of the control unit. The respective interface protocol can be obtained from the manufacturer.

11. Cleaning

Your unit should be cleaned regularly to ensure a long service life.

Cleaning of the unit is imperative in the following cases:

- if the piston unit has seized
- after titration of reagents whose use is restricted (see point 4 and 5)
- prior to changing reagent
- prior to sterilisation
- prior to replacing suction, ejection or recirculation valves
- prior to longer storage periods

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


Warning!

- Put on protective clothing, safety glasses and protective gloves prior to commencing cleaning.
 - Safety instructions (see point 1) must be observed.
 - The suction hose, ejection unit glass cylinder and valves are filled with liquid! All openings should be turned away from the body/from other persons.
-

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
11.1 Service/Cleaning

- Terminate titration, fit the closure cap and deactivate the unit.
- Place the unit mounted on a storage bottle into a cleaning container.
- Using suitable gloves, unscrew the unit from the storage bottle and remove (the suction hose no longer sticks into the liquid) **Warning: Liquid may emerge from the suction hose!**
- Hold the ejection unit over the storage bottle opening.
- Pull off the closure cap and push to the rear into the retention position.
- Activate the unit. Press “Venting” repeatedly until no more liquid emerges.

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
11.2 Basic cleaning with the venting function

- Immerse the suction hose in a suitable cleaning solution. Press the “Venting” button to purge the unit completely.


Caution:

- Residual liquid may cause splashing.
- Hold away from body at all times.

- Pull the suction hose out of the cleaning agent and purge the hose and unit with the “Venting” function until they are completely empty (until no more liquid emerges).
- Remove the following components from the unit:
 - Unscrew the ejector unit from the unit.
 - Pull off the recirculation hose.
 - Release the union nut on the suction hose and pull off the suction hose. Immerse these components in a suitable cleaning solution.

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11.3 Intensive cleaning

- Clean first as described in point 11.1. and 11.2.
- Remove the valve block, cylinder and piston.
- Turn the upper part of the unit up to the right-hand stop (the point is aligned with the right-hand edge of the inscription on the chrome ring), then draw 6 mm upwards and turn to the left up to the left-hand stop. Now pull the valve block and cylinder upwards and away from the upper part of the unit.

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
Note!

The piston can be removed by pushing it sideways away from the spindle if necessary. The spindle is held in position with a finger for this purpose, and the piston can be pushed to the side by hand to disconnect it from the spindle. For greater ease of access, the inspection window can also be removed if necessary.

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
11.3 Intensive cleaning

If the basic unit is in the cleaning position, position mounting can be facilitated by moving the spindle slowly to the lower stop by pressing the "Move" button in the message on the control unit for a longer period of time. If the spindle is already in the lowermost position, it can be raised easily by pressing the button again for a longer period of time.

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11.3 Intensive cleaning

- Push the upper part of the unit 6 mm onto the valve block and cylinder up to the first stop and turn to the right-hand stop, pushing a further 6 mm together subsequent to this.
- Now turn the ejection unit to the desired position.
- Check the function of the unit.

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Warning!

Insert the piston carefully into the glass cylinder, as the seal lip can otherwise be damaged. This will cause the system to leak!

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12. Sterilisation

The components in the unit conveying media are resistant to steam sterilisation up to 121 °C and an exposure time conforming to DIN $t_e = 20$ min.

Warning:

Only the following components should be steam sterilised:

- Valve block with cylinder
 - Suction hose
 - Ejection unit
 - Recirculation hose
 - Piston
-

The upper part of the unit, control unit and power supply unit should never be sterilised with steam!
(see point 2 for overview of individual components).

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12.1 Preparation

Cleaning should be realized prior to sterilisation (see point 11).

12.2 Dismantling




Note:

Dismantling can only be realized in the home position (the piston unit is in the lowest position). The unit should only be dismantled when deactivated!

In addition to the components listed in point 11.1. the other unit components conveying media described in point 11.2. are dismantled and cleaned.

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
12.3 Loosening the suction, ejection and recirculation valve

- Disconnect the suction hose.
- Disconnect the ejection unit.
- Loosen the suction valve with the assembly wrench by turning it counterclockwise a half a revolution.
- Loosen the ejection valve with the assembly wrench by turning it counterclockwise a half a revolution.
- Loosen the recirculation valve with the assembly wrench by turning it counterclockwise a half a revolution.

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
12.3 Loosening the suction, ejection and recirculation valve


- Grasp the cover laterally and pull off.

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- Open the nut in a counterclockwise direction by turning it half a revolution, and loosen the tap by applying pressure to the nut and dismantle the piston.

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
- Disconnect the piston from the spindle by pushing it to the side (see point 11.3).

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Important:

The piston should be cleaned in a suitable solution prior to sterilisation. Then push the piston into the glass cylinder to a point where the piston protrudes slightly out of the glass cylinder.

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
12.4 Sterilisation procedure

Realize steam sterilisation of all components conveying media (up to 121 °C, exposure time conforming to DIN te = 20 min).

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
Warning!

- Allow the unit to cool to room temperature and only then commence reassembly.
 - Check all components for damage and replace if necessary.
-

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12.5 Assembly

Remount the piston, valve block and cylinder:

Fit the upper part of the unit onto the valve and cylinder. Align the marking on the upper part of the unit with the left-hand edge of the marking on the chrome ring. Push the upper part of the unit 6 mm onto the valve block and cylinder up to the first stop, pushing a further 6 mm together subsequent to this. Now turn the ejection unit to the desired position.

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
Warning!


Insert the piston carefully in the glass cylinder, as the seal lip can otherwise be damaged. This will cause the system to leak!

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
12.6 Tightening the suction, ejection and recirculation valve


- Tighten the suction valve firmly with the assembly wrench by turning it in a clockwise direction.

- Tighten the ejection valve firmly with the assembly wrench by turning it in a clockwise direction.

- Tighten the recirculation valve firmly with the assembly wrench by turning it in a clockwise direction.

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
12.6 Tightening the suction, ejection and recirculation valve

- Tighten the nut on the valve tap by hand by turning it in a clockwise direction.

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12.7 Function inspection


- Activate the unit.
- Function inspection with the “Venting” symbol . The venting function should be realized (can even be realized without liquid).
- Please check “Trouble-shooting” in case of malfunction (see point 16) or contact the manufacturer.

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13. Valve replacement



Leaking or damaged suction, ejection and recirculation valves should be replaced immediately (see “Accessories and spare parts” in point 19)

13.1 Suction valve replacement

- Realize cleaning (see point 11).
- Unscrew the suction hose.
- Stand the unit on its head to prevent the valve ball falling out during replacement.
- Unscrew the suction valve with the assembly wrench and remove the valve ball.
- Insert the ball again and screw on the new suction valve, tightening with the assembly wrench.
- Push on the suction hose up to the stop and tighten the union nut firmly (see point 6.2).

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13.2 Ejection valve replacement

- Realize cleaning (see point 11).
- Unscrew the ejection unit.
- Screw out the ejection valve with the assembly wrench and remove.
- Screw in the new ejection valve and tighten firmly with the assembly wrench.
- Push on the ejection unit up to the stop and tighten the union nut firmly.
Check to ensure a secure fitting!

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13.3 Recirculation valve replacement

- Unscrew the recirculation valve with the assembly wrench and remove.
- Screw in the new recirculation valve and tighten firmly with the assembly wrench.

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14. Repair by manufacturer

- A defective unit should be cleaned carefully prior to sending for repair (see point 11).
- Units will only be repaired if they are cleaned by the user and are no longer contaminated with any substances which could pose a risk to health.
- Uncleaned units will be returned unrepaired to the sender!
- Describe the malfunction accurately.
- Specify the liquid used (important!).
- The guarantee is rendered null and void if the unit has been subject to intervention other than those described (e. g. application of force).

15. Volumetric inspection

The parameters specified in point 17 for the regularity (R%) and coefficient of variation (CV%) of a unit can be checked as follows:

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1. Weighing procedure (see for example DIN EN ISO 8655-6)

- Dispense the test volume with the unit (deionised water).
- Weigh the quantity dispensed on the analysis scale.
- Convert this quantity into volume.

Warning: consideration should be given to temperature and air pressure.

2. Repetitive weighing

The weighing procedure described in 1 should be repeated at least five times to obtain reliable details on the regularity and coefficient of variation of a unit.

3. Calculation the regularity and coefficient of variation

- The regularity (R%) and coefficient of variation (CV%) of the unit should be calculated on the basis of statistical quality control formulas:

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$$s = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n-1}}$$

$$CV\% = \frac{100 \cdot s}{\bar{V}} \times \frac{V_s}{V_o}$$

Conversion of the mass determined to volume

$$V = m \cdot Z$$

Mean value

$$\bar{V} = \frac{\sum_{i=1}^n V_i}{n}$$

Systematic measuring deviation (regularity R%)

$$R\% = \frac{100 \cdot \bar{V} - V_s}{V_o}$$

**Random measuring deviation
(coefficient of variation CV%)**

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Legend

m	=	weighed value determined in g
V	=	volume determined in ml
Z	=	correction factor (e.g.: 1.0029 µl/mg at 20 °C and 101.3 kPa 1.0040 µl/mg at 25 °C and 101.3 kPa)

\bar{V}	=	mean volume value
n	=	number of measurement
V_s	=	test volume
V_0	=	nominal volume
S	=	standard deviation

This brief guide merely illustrates the principal procedure and calculation for calibrating the unit. The precise procedure is described in detail in DIN EN ISO 8655.

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16. Trouble-shooting

16.1 Messages

The red flashing LED in the basic unit indicates that the unit is in a non-operational condition. The control unit also generates a suitable warning message:

Messages	Reason	Solution
Piston blockage, check basic unit	<ul style="list-style-type: none">- Recirculation lever not applied to stop- Valves are blocked- Liquid viscosity is too high- Closure cap is not removed	<ul style="list-style-type: none">- Turn the recirculation lever completely to the front or rear.- Clean or replace valves.- Select slower speed, use dispensing tip instead of titration tip.- Remove closure cap.
Unit is in cleaning position, valve block is not fitted	<ul style="list-style-type: none">- Upper part of basic unit is turned to the right-hand stop (dismantling position)- Upper part of basic unit has been dismantled	<ul style="list-style-type: none">- Turn the upper part slightly to the left.- Assembly basic unit again.

Messages	Reason	Solution
No basic unit identified	- Basic unit not connected to control unit	- Check cable connection.
No piston fitted	- Basic unit was assembled without piston	- Dismantle basic unit and fit piston again.

Special cleaning position case:

If the basic unit is in cleaning position, mounting can be facilitated by moving the spindle slowly to the lower stop by pressing the "Move" button in the message on the control unit for a longer period of time. If the spindle is already at the lowest point, pressing the button for a longer period will move it slightly upwards again.

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16.2 Other problems – solutions

If the unit can no longer be operated, it can be deactivated by pressing and holding (approx. 10 s) in the upper right-hand corner of the touch screen (battery symbol position) for a longer period. All changes realized since the last activation are lost as a result, so this function should only be used in emergencies.

Problem	Reason	Solution
Liquid emerges from the gap between the suction, ejection or recirculation valve and valve block	<ul style="list-style-type: none">- Valves loose or damaged	<ul style="list-style-type: none">- Tighten the valves firmly with the assembly wrench or replace if necessary. Warning: the suction, ejection, and recirculation valve should be tightened firmly with the assembly wrench after sterilisation.
Liquid emerges between the ejection valve and ejection unit	<ul style="list-style-type: none">- Ejection unit loose or damaged- Incorrect ejection unit	<ul style="list-style-type: none">- Mount original ejection unit correctly.- Use only original ejection unit.
Unit is suctioning in air	<ul style="list-style-type: none">- Unit not vented- Suction hose loose- Suction hose and/or union nut damaged	<ul style="list-style-type: none">- Vent unit.- Tighten suction hose firmly.- Replace suction hose.

Problem	Reason	Solution
Liquid emerges from recirculation valve (rear of unit)	<ul style="list-style-type: none"> - Recirculation valve loose or damaged 	<ul style="list-style-type: none"> - Tighten recirculation valve firmly with the assembly wrench or replace.
Liquid emerges between the valve shaft and valve block (only visible when cover is removed)	<ul style="list-style-type: none"> - Valve shaft loose 	<ul style="list-style-type: none"> - Tighten nut on valve shaft firmly.
Erroneous dispensing volume	<ul style="list-style-type: none"> - Ejection unit loose or damaged - Suction hose loose or damaged - Suction, ejection and recirculation valves loose or damaged 	<ul style="list-style-type: none"> - Mount original ejection unit correctly. Use only original ejection unit. - Tighten suction hose firmly. - Tighten the valves firmly with the assembly wrench or replace if necessary. Realize cleaning.

16.2 Other problems - solutions

Problem	Reason	Solution
Liquid over piston unit	- Damaged seal lip on piston unit	- Replace piston unit.
Visible air bubbles in glass cylinder	- Unit not vented	- Vent unit.
	- Suction hose or union not loose, deformed	- Tighten suction hose firmly. Change original suction hose.
Inaccurate dispensing volume	- Piston unit seal lip damaged	- Replace complete piston unit.
	- Ejection unit loose or damaged	- Mount ejection unit correctly.
Erroneous values on display	- Unit error	- Reset unit to default settings using the button in the "Main menu" item .

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16.3 Contacting the manufacturer

Should the measures as listed above fail to completely remedy the malfunction, the dealer or manufacturer should be contacted. Please specify the serial number and revision number of the basic unit and control unit. Serial numbers are embossed on unit components and can be viewed in the Info menu. Revision numbers can only be viewed in the Info menu.

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17. Technical data

Performance data Basic unit	10 ml	20 ml	50 ml
Unit type	Microprocessor-controlled, motorised piston lift system with recirculation		
Volume range	0.01 ml - 10.00 ml and 0.010 ml - 10.000 ml (max. piston lift 10 ml, automatic addition), resolution: 1 µl	0.1 ml - 20 ml (max. piston lift 20 ml, automatic addition), resolution: 10 µl	1.00 ml - 50 ml (max. piston lift 50 ml, automatic addition), resolution: 10 µl
Volume accuracy	R = 0.2%, CV = 0.07%	R = 0.2%, CV = 0.07%	R = 0.2%, CV = 0.05%
Velocity	in 16 stages from 0.01 ml/s to 2 ml/s independently adjustable for suctioning and dispensing speed	in 10 stages from 0.1 ml/s to 4 ml/s independently adjustable for suctioning and dispensing speed	in 10 stages from 1 ml/s to 10 ml/s independently adjustable for suctioning and dispensing speed
Operating temperature range	4 °C to 40 °C		
Power supply	Plug-in mains unit 15V, 1 A (optional via battery in control unit)		
Optional battery operation	Battery charge for an average of 2.5 to 3 h continuous operation. Battery capacity indicated on TFT display, charging time 6.5 h, intelligent electronics prevent overloading and deep discharging.		
Materials	Housing made of PA6 reinforced glass fibre, all parts in contact with liquid are made of chemically inert material		
Cleaning, sterilisation	Motorised rinsing function, can be disconnected without tools, media-conveying components are suitable for steam sterilisation (121 °C)		
Aeration	Automatic at the push of a button		
Dimensions/Weight	320 mm x 80 mm x 63 mm (H x D x W)/ 660 g		
Quality control	in compliance with DIN EN ISO 8655		
Documentation	DIN EN ISO 9001 compliant certification through unit protocol with serial number		

17. Technical data

Performance data Control unit	10 ml	20 ml	50 ml
Module type	External control module incl. touch screen control		
Touchscreen	3.5" with 320 X 240 pixel, most important functions directly accessible from main screen		
Titration volume preselection	Up to 10 ml	Up to 20 ml	Up to 50 ml
Rapid titration	For rapid medium dispensing at beginning of titration process		
Droplet titration	With 1 droplet per button push, accelerating if button is pressed		
Quick-Cal	Calibration of other temperature or medium possible		
Method memory	Up to 9 different unit configurations can be saved		
Languages	4 languages (D, GB, F, S) integrated as standard option		
PC interface	USB or RS232		
Remote control	Unit can be remotely controlled via an additional cable-connected device		
Dimensions/Weight:	55 mm x 130 mm x 110 mm (H x D x W) /230 g		
Scope of delivery			
Unit with A 45 internal thread and additional A 32, A 38 and S 40 threaded adapters			
Suction hose with union nut, recirculation hose, titration ejection unit, valve fitting wrench, touchscreen control unit "titration", brown light protection window, power supply unit, data power cable, operating manual and individual quality certificate			

The power supply unit may deviate from the specific national design in the scope of delivery.

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Conformity

Conformity certificate conforming to German calibration regulations (appendix 12).

CE

This symbol confirms that the unit meets the stipulated requirements of pertinent EU directives and is tested in conformance with these.

18. Materials

Piston	PTFE Compound
Cylinder	Borosilicate Glass 3.3
Valve block	Halar (ECTFE)
Valve mantle	Halar (ECTFE)
Valve ball	High-purity aluminium oxide ceramic (99.7%)
Valve seat	High-purity aluminium oxide ceramic (99.7%)
Valve spring	Hastelloy

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for unit(Serial number)

Medium used:

	YES	NO	Please define in case of YES
Blood, bodily fluids, patholog. samples	<input type="checkbox"/>	<input type="checkbox"/>	_____
Other hazardous biologic substances	<input type="checkbox"/>	<input type="checkbox"/>	_____
Biologically-degradable materials which can pose a risk	<input type="checkbox"/>	<input type="checkbox"/>	_____
Chemicals containing substances harmful to health	<input type="checkbox"/>	<input type="checkbox"/>	_____
Radioactive materials	<input type="checkbox"/>	<input type="checkbox"/>	_____
Other hazardous and toxic substances	<input type="checkbox"/>	<input type="checkbox"/>	_____

We hereby declare that

- the unit has been carefully cleaned and decontaminated prior to shipment,
- no risks resulting from contamination exist,
- we are aware that liability claims can be made against us in the event of damage being caused by contaminated units.

Address of sender:

Name: _____

Date: _____

Signature: _____