# HIRSCHMANN®

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## Instruction manual

Please read before use!

# rotarus®

Handhabung Einkanal Pumpenkopf  
singlechannel pump head handling

Hirschmann Laborgeräte GmbH & Co. KG  
Hauptstraße 7-15 • 74246 Eberstadt Germany  
Fon +49 7134 511-0 • Fax +49 7134 511-990  
[www.hirschmannlab.com](http://www.hirschmannlab.com) • [info@hirschmannlab.com](mailto:info@hirschmannlab.com)

X2120437000

## Inhalt

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## Content

	Page
1. Pump head assembly .....	3
2. Selection and insertion of the pump tube.....	6
2.1 Selection of tube material.....	7
2.2 Tube insertion .....	7
2.3 Tube fitting for pump drives with automatic tube identification .....	9
2.4 Tube care and exploitation of maximum service life .....	10
2.5 Running-in phase for tubes.....	11
2.6 Pumping against pressure .....	11
2.7 Flow rates.....	12

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
## 1. Pump head assembly

The standard, flow and volume unit versions have a pump head quick-change system.

All single channel pump heads are fitted to the drive in the same manner.

Pump heads from other manufacturers can be used with appropriate adapter plates. Fitting is realized in the same manner.

The cover on the pump drive should be removed prior to mounting the pump head. Pull the pin on the lower part of the pump mounting ring downwards for this purpose, then remove the cover with the other hand. See ill. p. 5.

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## 1. Pump head assembly

The cover removed from the pump should be fitted to the pump head. For reasons of safety, the single channel pump head should not be operated without this cover, as there is otherwise a risk of injury at the pump head adapter position when the pump is running.

Fit the pump head to the pump drive and push against the pump until it engages.


### **Assembly made easier**

- For flow and volume:  
Actuating the MIN/MAX button causes the motor to run at low speed until the pump drive and pump head coupling components are congruent.
- For standard:  
a lower speed should be selected here.

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## Pump head dismantling/changing

Stop or deactivate the pump.

Pull the pin on the lower part of the pump mounting ring downwards and pull the pump head forwards with the other hand.

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---

## 2. Selection and insertion of the pump tube

---

### **Warning!**

Only tubes may be used which have the correct dimensions and characteristics suited to the pump heads used. Significant parameters for reliable function are the wall thickness and tube material hardness. Incorrect tubes can negatively impair the function or lead to damage to the tube or pump.

Only tubes with a wall thickness of 1.6 mm and an inner diameter of between 1.6 mm and 11.2 mm should be used for the PK 10-16 single channel pump head.

Only tubes with a wall thickness of 2.4 mm and an inner diameter of between 4.8 mm and 9.6 mm should be used for the PK 10-24 single channel pump head.

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---

## 2.1 Selection of tube material

Ensure that the tube material to be used is resistant to the medium being conveyed.

The user is responsible for ensuring use of the correct tube material.

## 2.2 Tube insertion

---

### **Warning!**

Liquid may be sprayed out by the pump in the event of the tube rupturing. Appropriate measures should be taken to protect operating personnel and equipment.

Switch off the drive prior to inserting the tube.

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
**Warning!**

Incorrect combinations lead to malfunctions and premature wear.

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
Press the small lever on the right-hand side forwards.

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
Close the tubing bed and press downwards until it engages (this is the precondition for starting the pump in the flow and volume variants).

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

### 2.3 Tube insertion for pump drives with automatic tube identification

In the case of pump tubes with a tag:  
Insert the tag of one of the tubes used into the slot provided on the pump head.

In the case of pump tubes without a tag or pump heads from other manufacturers, data must be programmed manually using the control unit (see instruction manual for tube and head monitoring).

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## 2.4 Tube care and exploitation of maximum service life

We recommend relieving tension on the tube or removing it from the pump head in the event of longer operating interruptions.

Opening of the tubing bed is adequate in the case of a single channel pump head.

### **Warning!**

The medium can flow back where the tube is relaxed.

The tube service life depends on the tube material and application conditions and is governed by the chemicals used, selected speed, differential pressure, temperature, viscosity, etc. Non-binding reference values can be found in the tube descriptions.

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## 2.5 Running-in phase for tubes

Every new tube changes its behaviour during the first operating hours. It requires a certain running-in period to adjust.

It is necessary to condition the new tube for 1-3 hours prior to use with water or the medium to be pumped to achieve a consistent and reproducible flow rate.

## 2.6 Pumping against pressure

A pressure of approx. 1.5 bar can be achieved with the peristaltic pump. Slightly higher pressure can also be achieved for brief periods, depending on the pump head, tube material and tube dimensions used.

Please contact the pump manufacturer in cases of doubt.

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## 2.7 Flow rates

See data sheet for respective pump head.

This can be found at  
[www.hirschmannlab.com/en/Info-Center/Downloads/Datenblaetter](http://www.hirschmannlab.com/en/Info-Center/Downloads/Datenblaetter)