HIRSCHMANN®

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Instruction manual

Please read before use!

rotarus®

Handhabung Mehrkanal Pumpenkopf
multichannel pump head handling

Hirschmann Laborgeräte GmbH & Co. KG
Hauptstraße 7-15 • 74246 Eberstadt Germany
Fon +49 7134 511-0 • Fax +49 7134 511-990
www.hirschmannlab.com • info@hirschmannlab.com

X2120438000

Inhalt

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Content

	Page
1. Pump head assembly	3
2. Selection and insertion of the pump tube.....	5
2.1 Selection of tube material.....	6
2.2 Tube insertion	6
2.3 Tube and cassette fitting.....	8
2.4 Tube and cassette fitting for pump drives with automatic tube identification	12
2.5 Removing the tube cassettes.....	13
2.6 Tube care and exploitation of maximum service life	14
2.7 Running-in phase for tubes.....	15
2.8 Pumping against pressure	15
2.9 Flow rates.....	16

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1. Pump head assembly


The standard, flow and volume unit versions have a pump head quick-change system.

All multichannel pump heads are fitted to the drive in the same manner.

Pump heads from other manufacturers can be used with appropriate adapter plates. Fitting is realized in the same manner.

The cover on the pump drive should be removed prior to mounting the pump head. Pull the pin on the lower part of the pump mounting ring downwards for this purpose, then remove the cover with the other hand.

Fit the pump head to the pump drive and push against the pump until it engages.

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
Assembly made easier

- For flow and volume:
Actuating the MIN/MAX button causes the motor to run at low speed until the pump drive and pump head coupling components are congruent.
- For standard:
a lower speed should be selected here.

Pump head dismantling/changing

Stop or deactivate the pump.

Pull the pin on the lower part of the pump mounting ring downwards and pull the pump head forwards with the other hand.

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

2. Selection and insertion of the pump tube

Warning!

Only tubes may be used which have the correct dimensions and characteristics suited to the pump heads used. Significant parameters for reliable function are the wall thickness and tube material hardness. Incorrect tubes can negatively impair the function or lead to damage to the tube or pump.

Only 2-stopper tubes (with a distance of 104 mm between stoppers) can be used for the cassettes.

Use of tubes without stoppers leads to malfunctions.

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2.1 Selection of tube material

Ensure that the tube material to be used is resistant to the medium being conveyed.

The user is responsible for ensuring use of the correct tube material.

2.2 Tube insertion

Warning!

Liquid may be sprayed out by the pump in the event of the tube rupturing. Appropriate measures should be taken to protect operating personnel and equipment.

Deactivate the drive before removing or inserting tube cassettes. The safety features provided should prevent fingers being inserted between the roller mechanisms and the base of the module. Nevertheless, care should be exercised in the interest of safety and fingers kept away from this area.

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2.2 Tube insertion

Warning!

Use only appropriate cassettes and suitable stopper tubes for the pump heads. Incorrect combinations lead to malfunctions and premature wear.

The **small cassette** is only suitable for use with the MKF 4-8 and MKF 12-8 pump heads.


The **medium cassette** is only suitable for use with the MKF 8-4 pump head.

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
2.3 Tube and cassette fitting

1. The tubing bed should be set to the maximum opening prior to inserting the tube in the cassette. This is achieved by turning the adjusting screw in an anticlockwise direction.

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2. Insert the stopper tube into the cassette so that the stoppers on both ends are resting on the tube bracket. The stopper should make contact with the tube holder.

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

2.3 Tube and cassette fitting

3. Suspend the cassette fitted with the pump tube in the first pump head carrier rod while holding the tube, fold over the roller head and engage in the second carrier rod.

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4. Optimum adjustment of contact pressure



Increase the contact pressure while the pump is running by turning the adjusting screw in a clockwise direction until the maximum flow rate is achieved. Then reduce the contact pressure slightly by turning the adjusting screw in an anticlockwise direction without changing the flow rate. Use of a lower speed facilitates adjustment. This procedure enables you to achieve the maximum pump tube service life.

5. Where the pump head is only partially fitted with cassettes, the foremost cassette can be fixed in position with the O-rings on the two carrier rods.

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
2.4 Pump drives with automatic tube identification

In the case of pump tubes with a tag:
Insert the tag of one of the tubes used into the slot provided
on the pump head.

In the case of pump tubes without a tag
or pump heads from other manufacturers,
data must be programmed manually
using the control unit (see instruction
manual for tube and head monitoring).

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2.5 Removing the tube cassettes

1. Stop the pump.
2. Press lightly from above on the cassette on the engaging side, pull the engaging lever upwards and pivot the cassette upwards.
3. Remove the cassette.

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2.6 Tube care and exploitation of maximum service life

We recommend relieving tension on the tube or removing it from the pump head in the event of longer operating interruptions.

Releasing the cassettes on one side is adequate in the case of multichannel pump heads.

Warning!

The medium can flow back where the tube is relaxed.

The tube service life depends on the tube material and application conditions and is governed by the chemicals used, selected speed, differential pressure, temperature, viscosity, etc. Non-binding reference values can be found in the tube descriptions.

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2.7 Running-in phase for tubes

Every new tube changes its behaviour during the first operating hours. It requires a certain running-in period to adjust.

It is necessary to condition the new tube for 1-3 hours prior to use with water or the medium to be pumped to achieve a consistent and reproducible flow rate.

2.8 Pumping against pressure

A pressure of approx. 1.5 bar can be achieved with the peristaltic pump. Slightly higher pressure can also be achieved for brief periods, depending on the pump head, tube material and tube dimensions used.

Please contact the pump manufacturer in cases of doubt.

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2.9 Flow rates

See data sheet for respective pump head.

This can be found at
www.hirschmannlab.com/en/Info-Center/Downloads/Datenblaetter