# HIRSCHMANN®

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## Instruction manual

Please read before use!

# solarus®

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## Content

	page
1. Warnings and safety instructions .....	4
2. Individual components .....	8
3. Principle of operation .....	9
4. Prohibited applications .....	10
5. Restricted applications .....	11
6. Preparation of device .....	12
6.1 Using the light protection window .....	12
6.2 How to connect the suction tube .....	13
6.3 How to connect the recirculation tube .....	14
6.4 How to connect the discharge tube .....	14
6.5 Secure device on bottle .....	16
6.6 Align device with bottle label .....	16
7. Titrating .....	17
7.1 Remove closure cap .....	17
7.2 How to prime the device .....	18
7.3 How to fill the digital burette .....	19
7.4 Titration .....	20
7.5 Display modification .....	21
7.6 Customer-specific calibration Quick-Cal .....	22
8. Cleaning .....	24
8.1 Maintenance / cleaning .....	25
8.2 Thorough cleaning .....	29
9. Sterilization .....	32
9.1 Preparations for sterilization .....	32

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	page
9.2 Disassembly .....	33
9.3 How to loosen the suction, discharge and recirculation valve .....	33
9.4 How to perform a sterilization .....	35
9.5 How to reassemble the piston .....	35
9.6 How to secure suction, discharge and recirculation valve .....	37
9.7 How to assemble the device .....	38
10. How to exchange valves .....	39
10.1 How to exchange the suction valve .....	39
10.2 How to exchange the discharge valve .....	41
10.3 How to exchange the recirculation valve .....	42
11. Repairs at the manufacturer .....	42
12. How to test the correct volume .....	44
13. Trouble-shooting .....	48
14. Technical data / materials .....	50
15. Accessories - spare parts .....	53
16. Clearance certificate form .....	56

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---

### **Please note:**

Please read the user manual and safety instructions carefully. Each person working in the laboratory must be familiar with these safety instructions prior to using the device. The manual must be accessible at any time. Not every possible safety issue is covered in this manual. Each applicant of the device is solely responsible to maintain safety and health standards and to determine and follow restrictions prior to using the device.

- 
- All general safety instructions of the laboratory such as wearing protective clothing, safety goggles and gloves appropriate for the application of certain fluids must be followed at all times.

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## 1. Warnings and safety instructions

- Instructions provided by the reagent manufacturer must be observed at any time.
- This device is only intended for titrating fluids. Certain applications are either prohibited or restricted. These instructions must be strictly adhered to (If in doubt contact the manufacturer).
- Prior to use, all connections must be checked for leakage. The correct positioning of the discharge spout and the suction tube must be verified. If the discharge spout is not secured properly fluid could leak or splatter out during the operation of the burette.
- The closure cap must be removed from the discharge tube prior to titrating. Hand wheels should not be turned as long as the closure cap is still positioned on the device.

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## 1. Warnings and safety instructions


- Make sure that you do not endanger yourself or any other person. Never work towards another person. Avoid splattering. Use the appropriate receiving vessels. Always titrate onto the inner wall of the receiving container instead of its content or onto the bottom of the container.
- Never use force!
- As soon as a titration has been completed position closure cap on discharge tube. Caution! Reagent may leak out!
- Use original accessories and replacement parts only.
- Once sterilization has been completed (refer to section 9) let the digital burette cool down to room temperature before use.

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## 1. Warnings and safety instructions

- In the event of errors, such as
  - slow moving or jammed piston
  - sticky or leaking suction/discharge/recirculation valve
  - leakage
  - broken parts
  - Tubes
  - Sealing ringsstop working immediately. Before you continue to use the device read the instructions for cleaning of the device (section 8) and trouble-shooting (section 13). If necessary send the unit to the manufacturers' service center for repair. Please read section 11 for more information.
- Volume dispensing complying with standards can only be guaranteed if the standard accessories provided are used.
- Where the instrument screwed to the bottle is to be carried to another working location, refrain from lifting the combination by holding the instrument alone. The bottle should also be securely held during carrying.

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## 2. Individual components

1. Discharge tube complete
2. Recirculation valve
3. Suction valve
4. Discharge valve
5. Stopcock with recirculation valve control,  
nut, washer, O-ring
6. Cover

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### 3. Principle of operation

- The device is a manually operated bottle top burette with digital display for the titration of fluids. The basic unit has thread A45 built in. Multiple adapters (refer to accessories) permit the use of standard bottles.
- Hold the hand wheels with your thumbs and index fingers, move your thumbs upward, turning the hand wheels away from the body, the suction valve opens and fluid is drawn into the cylinder of the digital burette via the suction tube from the bottle. As soon as the suction process is either interrupted or completed the suction valve closes.

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### 3. Principle of operation

- Once the burette has been filled, hold the hand wheels with the thumbs and index fingers, move your thumbs downward, turning the hand wheels towards the body, the discharge valve opens and dispenses fluid from the cylinder through the discharge tube. The titrated volume is shown in the LC-Display of the digital burette as the process proceeds.

### 4. Prohibited applications

---

#### **The digital burette may not be used with:**

- Hydrofluoric acid, since it attacks glass
  - Fluids that attack glass, Halar (ECTFE), FEP and Hastelloy
  - Suspensions, since solid particles may clog or block valves
  - Easily-combustible media (note the flashpoint)
  - Use in a potentially-explosive environment
-

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## 5. Restricted applications

- highly concentrated acids, alkalines, salt solutions, HF.
- the burette is suitable for operating temperatures between +4°C and +50°C.
- The unit should only be used in areas in which it can be assured that no hazardous potentially-explosive atmosphere will occur.
- Do not deposit in potentially-explosive areas.

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## 6. Preparation of the device

---

### **Attention, please note prior to assembly!**

- Proper and safe operation is only warranted if the enclosed original parts are used.
  - Never use aftermarket or damaged parts.
- 

### 6.1 Using the light protection window



The existing clear inspection window can be replaced with the colored light protection window to protect media which are sensitive to light.

The clear inspection window should be raised initially on one side for removal to ensure that it slides out of its fixture, and then raised again on the other side.

The colored light protection window can now be suitably fitted (note correct position).

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
## 6.2 How to connect the suction tube

- Hold suction tube next to supply bottle with the retainer ring pointing towards the bottom of the bottle.
- Cut the tube at an angle at the height of the lower edge of the rim on the reservoir bottle.
- Attach suction tube securely.
- Next screw retainer ring onto suction valve.

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### 6.3 How to connect the recirculation tube


- Firmly press recirculation tube in the designated recirculation outlet of valve block.

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### 6.4 How to connect the discharge tube


- Discharge tube includes:
  - tube guide
  - retainer ring
  - dispensing tube
  - closure cap

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
#### 6.4 How to connect the discharge tube

- Attach discharge tube securely.

- Firmly tighten retainer ring and verify correct positioning.

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
#### **Please note:**

- Confirm correct positioning of the discharge tube.
- Only use original discharge tube (refer to accessories and replacement parts).
- Never install a deformed or defective discharge tube.

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### 6.5 Secure device on bottle

- Screw device either with or without adapter onto reservoir bottle.

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### 6.6 Align front of device with bottle label

- Rotate device to desired position of bottle label

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---

#### Warning!

- Always follow all safety instructions and avoid prohibited and restricted applications.
  - Always wear protective clothing, goggles and gloves.
-

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
---

### Please note:

- Always follow safety instructions (refer to section 1) as well as application prohibitions and restrictions (refer to sections 4 and 5).
  - Only use digital burette in such a way that the discharge tube is never directed towards the user or other people.
  - Do not press piston down when closure cap is in place.
- 

### 7.1 Remove closure cap


- Place receiving container underneath discharge tube.
- Remove closure cap as shown in the figure.

**Please note: Reagent may drip during this process!**

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
### 7.1 Remove closure cap

- Press closure cap back on as shown.
- During operation the closure cap should be placed in the safe storage position.

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
### 7.2 How to prime the device

During the initial filling of the digital burette or in case air bubbles should occur in the cylinder chamber, the recirculation valve control must be set to point towards the back of the device. Then hold the hand wheels with your thumbs and index fingers, move the thumbs downward, turning the hand wheels towards the body until the piston is moved down completely. Turn the hand wheels a few times briefly in both directions.

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


## 7.2 How to prime the device

Make sure that the piston is moved completely to the bottom. When bubbles are no longer visible in the view glass the recirculation valve control should be set to titration (pointing toward the discharge tube) and the actual titration may begin.

## 7.3 How to fill the digital burette

Touch one of the two keys briefly to switch on the device. Press the left key to switch the device to the „Fill“ mode (arrow points up in the display). The volume indication in the display is then automatically set to zero. - Fill the burette by holding the hand wheels with your thumbs and index fingers, move your thumbs upwards, turning the hand wheels away from the body.

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The burette may either be filled completely (10 ml, 20 ml or 50 ml) or only partially.

## 7.4 Titration

To switch on the device touch one of the two keys briefly. Press the left key to switch the device to the „Titration“ mode (arrow points down in the display). The last volume dispensed appears in the display. The displayed value can be reset to zero by pressing the key on the right. Holding the hand wheels with the thumbs and index fingers, moving your thumbs downwards, turning the hand wheels towards the body starts the dispensing of the titrant. The manufacturer suggests

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#### 7.4 Titration

dispensing the first few drops of the titrant into a waste container prior to the actual titration. To begin the titration process reset the display value to zero. Once the titration has been completed the displayed value is to be recorded.

#### 7.5 Display modification

The volume indicator is not reduced by turning up the hand wheel, but remains at its current value. While turning up the hand wheel, the upper “Mode arrow” flashes. To activate and de-activate this feature, press and hold the left key for 6 seconds.

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## 7.6 Customer-specific calibration Quick-Cal

The default calibration for the digital burette is done at 20 °C according to regulation DIN EN ISO 8655. If necessary the device may easily be calibrated to customer-specific applications (other temperatures etc.). To do this the right key must be pressed for 4 sec. The word „cal“ will then be flashing in the display as well as the current calibration factor.

The default setting is 10.00 for a 10 ml device, 20.00 for a 20 ml device and 50.00 for a 50 ml device. The calibration factors may be set between 9.00 and 10.99 for a 10 ml device, between 18.00 and 21.99 for a 20 ml device and between 45.00 and 54.99 for a 50 ml device. Pressing the right key modifies the actual decimal position flashing at the time. When the left key is pressed the device switches to the next decimal position.

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## 7.6 Customer-specific calibration Quick-Cal

After calculating the actual volume of the dispensed fluid (ideally use the maximum burette volume), the new factor can be entered directly as the customerspecific calibration factor (higher factor means smaller volume output). The new factor should be determined using the actual weight of the delivered fluid and the specific conditions in the user laboratory (temperature, atmospheric pressure etc.). This should be done in accordance with GLP.

Press the left key for 4 sec and the selected value will be stored. The display „Pro“ will appear for a moment. On the left margin of the display the word „cal“ indicates, that the device is operating with a customer-specific calibration.

To switch between the customer-specific calibration („cal“ appears in display) and the default calibration („cal“ does not appear in the display) press the left key for 4 sec.

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## 8. Cleaning

To warrant the lifetime of your digital burette, the device must be cleaned on a regular basis.

The following situations require the device to be cleaned:

- After titration of reagents whose application is restricted (refer to section 5)
- Prior to change of reagent
- Prior to sterilization
- Prior to exchange of suction, discharge and recirculation valves
- Prior to longer storage periods

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
### **Please note:**

- Before the attempt is made to clean the device, put on protective clothing, goggles and gloves.
- Safety instructions (refer to section 1) must be observed at all times.
- Suction tube, discharge tube, cylinder and valves are filled with fluids. Direct all openings away from your body and other persons.

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### **8.1 Maintenance / cleaning**


- Complete titration and put closure cap in place (refer to section 7.5).
- Place the device, which is attached to a reservoir bottle into the cleaning vessel.

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
- Unscrew and remove device from reservoir bottle wearing suitable hand protection. (Suction tube should be out of the fluid in the bottle).

### Please note:

Fluid may spill out of the suction tube.

- Hold discharge tube above the opening of the reservoir bottle.
- Remove closure cap and park in storage position.
- Either discharge all of the remaining fluid from the digital burette back into the supply bottle, or dispose it in accordance with the GLP.

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
## 8.1 Maintenance / cleaning

- Completely empty the device by rotating the hand wheel back and forth several times.
- Dip suction tube into a suitable cleaning solution and rinse the device thoroughly by drawing and discharging the cleaning solution several times.


### **Caution:**

Remaining fluids may splatter!  
Always direct the spout away from your body!

- Rinse digital burette with cleaning solution (e.g. distilled water or acetone) through repeated filling and discharging of fluid.
- Remove suction tube from solvent and empty device completely through repeated piston movement.

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
- Remove the following parts from the burette as described below and place in cleaning solution:
- Unscrew discharge tube from digital burette

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

## 8.2 Thorough cleaning

- First perform cleaning as described in section 8.1 above.
- Then remove valve block, cylinder and piston:
- Rotate upper section of the device to the right until the movement stops (arrow is lined up with the right edge of the chrome ring), then pull out upper section 6 mm and turn to the left until movement stops. Remove valve block with cylinder from upper section by pulling downwards.
- If necessary, completely remove piston and toothed rack from the upper section of the device by rotating the hand wheels.

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
**Please note!**

If necessary, piston and toothed rack can be separated by sliding the piston of the rack.

- Check piston for damage on the sealing lips (exchange if needed), then place piston into cleaning solution as described in section 8.1.
- Clean cylinder with brush and cleaning solution, then rinse with rinsing solution.
- Reassemble piston, toothed rack, valve block and cylinder.
- Slide piston onto toothed rack.  
For reassembly, align the toothed rack in such a way, that the teeth face the view glass. Then turn the hand wheels in the fill direction and guide the rack into the gears.

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- Place upper section of the device onto the valve block with cylinder. The arrow on the upper part must align with the left edge of the marking on the chrome ring. Insert upper section 6 mm into valve block with cylinder until top, then rotate to the right until movement stops and push down another 6 mm. Move discharge spout to the desired working position.

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### Please note:

- Carefully insert piston into the glass cylinder so that sealing lip will not be damaged. A damaged sealing lip will cause leakage problems!
- Verify proper functioning of the device.

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## 9. Sterilization

All parts of the digital burette which come in contact with fluid can be steam-sterilized up to 121 °C/2 bar and sustain a sterilization time according to DIN te = 20 min.

### 9.1 Preparations for sterilization

- The burette must be cleaned (refer to section 8.1) prior to sterilization.

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## 9.2 Disassembly


- In addition to the parts listed under 8.1 all other parts of the burette which come in contact with fluid must be disassembled and cleaned as described under 8.2.

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## 9.3 How to loosen suction, discharge and recirculation valve


- Loosen suction valve by turning it 180 degrees counter-clockwise with enclosed wrench.

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

### 9.3 How to loosen suction, discharge and recirculation valve


- Loosen discharge valve by turning it 180 degrees counter-clockwise with enclosed wrench.

- Loosen recirculation valve by turning it 180 degrees counter-clockwise with enclosed wrench.

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
## 9.4 How to perform a sterilization

Perform steam-sterilization of all parts, coming in contact with fluid (up to 121 °C/ 2 bar, sterilization time according to DIN te = 20 min.).

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### Please note:

- Allow device to cool down to room temperature, then begin reassembly.
- Check all parts for damage and replace if needed.

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
## 9.5 How to reassemble the piston

- Reassemble piston, toothed rack, valve block and cylinder:  
Slide piston onto toothed rack.

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For assembly, align the toothed rack in such a way, that the teeth face the view glass. Then turn the hand wheels into the fill direction and guide the rack into the gears.

- Place upper section of the device onto the valve block with cylinder. The arrow on the upper part must align with the left edge of the marking on the chrome ring. Slide upper section 6 mm onto valve block with cylinder until stop, then rotate to the right until movement stops and push down another 6 mm. Finally move discharge tube to the desired working position (arrow should no longer be aligned with the right edge of the chrome ring) to lock the mechanism.

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
**Please note:**

Carefully insert piston into the glass cylinder so that sealing lips will not be damaged. A damaged sealing lip will cause leakage problems!


- Verify proper functioning of the device.

### 9.6 How to secure suction, discharge and recirculation valve

- Tighten suction valve clockwise with wrench.
- Tighten discharge valve clockwise with wrench.
- Tighten recirculation valve clockwise with wrench.

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#### 9.6 How to secure suction, discharge and recirculation valve

- Secure nut on valve shaft clockwise by hand.
- Position cap.

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#### 9.7 How to assemble the device

- Assemble device (refer to section 6) and verify proper functioning (refer to section 7). In the event of malfunctioning refer to trouble-shooting (section 13).

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## 10. How to exchange valves


Leaking or damaged suction, discharge and recirculation on valves must be exchanged immediately (refer to accessories and replacement parts in section 15).

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### 10.1 How to exchange the suction valve


- Clean burette (refer to section 8).
- Disconnect suction tube.
- Turn device upside down. Otherwise valve ball may fall out during the exchange.
- Remove suction valve including valve ball with the help of the wrench.

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#### 10.1 How to exchange the suction valve

- Insert ball and install new suction valve. Tighten with wrench.

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

- Reattach suction tube snugly. Next screw retainer ring onto suction valve (refer to section 6.2).

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## 10.2 How to exchange discharge valve

- Clean burette (refer to section 8).
- Disconnect discharge tube (refer to section 8).
- Unscrew discharge valve with wrench and remove.
- Screw in new discharge valve and tighten with wrench.
- Position discharge tube tightly and secure retainer ring (refer to section 6.4). Check for proper positioning!

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### 10.3 How to exchange the recirculation valve

- Unscrew recirculation valve with wrench and remove.
- Screw in new recirculation valve and tighten with wrench.

### 11. Repairs at the manufacturer

- Prior to shipping a defective device to the manufacturer for repair the unit must be carefully cleaned (refer to section 8).
- Devices will only be repaired if they have been properly cleaned and decontaminated.

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#### 11. Repairs at the manufacturer

- Without proper proof of decontamination, devices will be considered contaminated. They may be serviced at higher costs or returned unserviced.
- The malfunctioning of the device should be described in detail.
- Please indicate which reagents have been used (Important!).
- If the device has been applied other than for its intended use or if the use of force is evident any claims under warranty are void.

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## 12. How to test the correct volume

The parameters listed under section 14 for accuracy (R%) and coefficient of variation (CV%) of a device can be tested as follows:

### 1. Weighing process (refer to DIN 8655-6):

---

- Fill the burette completely with de-ionized water and discharge the nominal volume.
- Weigh discharge quantity on analytic scale.
- Convert quantity into volume.

---

#### **Please Note:**

Temperature, air pressure and material must be taken into consideration.

---

### 2. Repeated weighing process

---

- The weighing process described above should be repeated at least five times to obtain reliable test results for the accuracy as well as coefficient of variation of a device.

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### 3. Computation of accuracy and coefficient of variation

- Accuracy (R%) and the coefficient of variation (CV%) of the device must be computed according to the formulas of statistical quality control.

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### 13. Trouble-shooting

Symptom	Possible cause	Solution
Fluid leaks at the gap between the valve block and the valves	<ul style="list-style-type: none"><li>- Valve loose or damaged</li></ul>	<ul style="list-style-type: none"><li>- Tighten valves with the wrench or replace if necessary (see section 10)</li></ul> <p><b>Note:</b> After sterilization (see section 9), suction/discharge/recirculation valves need to be tightened.</p>
Fluid is leaking between the discharge tube and the valve	<ul style="list-style-type: none"><li>- Discharge tube is loose or damaged</li><li>- Wrong discharge tube used</li></ul>	<ul style="list-style-type: none"><li>- Mount the original discharge tube correctly (see section 6.4).</li><li>- Use only the original discharge tube.</li></ul>
Device aspirates air	<ul style="list-style-type: none"><li>- Device not primed</li><li>- Suction tube loose</li><li>- Suction tube or retainer ring damaged</li></ul>	<ul style="list-style-type: none"><li>- Prime the device (see section 7.2).</li><li>- Tighten suction tube (see section 6.2).</li><li>- Replace section tube (see accessories spare parts section 15).</li></ul>

### 13. Trouble-shooting

Fluid leaks at the gap between the valve block and the recirculation valve	<ul style="list-style-type: none"><li>- Recirculation valve loose or damaged</li></ul>	<ul style="list-style-type: none"><li>- Tighten recirculation valve with the wrench or replace if necessary (see section 10)</li></ul>
Fluid leaks between stopcock shaft and valve block (only perceptible when cover is removed)	<ul style="list-style-type: none"><li>- Stopcock shaft loose</li></ul>	<ul style="list-style-type: none"><li>- Tighten nut on stopcock shaft</li></ul>
Delivered volume incorrect	<ul style="list-style-type: none"><li>- Discharge tube loose or damaged</li><li>- Suction tube loose or damaged</li><li>- Suction/ discharge/ recirculation valves loose or damaged</li></ul>	<ul style="list-style-type: none"><li>- Mount original discharge tube correctly (see section 6.4)</li><li>- Tighten suction tube (see section 6.2)</li><li>- Tighten valves with the wrench</li><li>- Clean the device (see section 8.1)</li></ul>
No display	<ul style="list-style-type: none"><li>- Device switched off</li><li>- Inadequate light / Solar cell concealed</li></ul>	<ul style="list-style-type: none"><li>- Briefly touch either one of the keys to switch on the device</li><li>- Find a brighter light source</li></ul>

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Digitalburette Volumen ml	R* %	CV* %
10 ml	* 0,3	* 0,1
20 ml	* 0,2	* 0,1
50 ml	* 0,2	* 0,1

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\* Indicated values for accuracy (R) and coefficient of variation (CV) relative to the maximum volume at a

- consistent temperature (20 °C) of the device, ambient air and fluids (de-ionized water)
- even and jerk-free operation

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#### 14. Technical data / materials

The digital burette has a test certificate, which certifies the conformity with the German weights and measures regulations (attachment 12).

Ambient temperature range for proper operation of the device and fluids: +4 °C and +50 °C.

**Material****Materials**

Kolben	Zylinder	Ventil- block	Ventil- gehäuse	Ventil- kugel	Ventil- sitz	Ventil- feder	Schläuche	Dichtungs- ringe
piston	cylinder	valve block	valve mantle	valve ball	valve seat	valve spring	tubes	sealing rings

**Digital-  
bürette**

PTFE Compound Halar (ECTFE)	Borosilikat- Glas 3.3	Halar (ECTFE)	Halar (ECTFE)	hochreine Al <sub>2</sub> O <sub>3</sub> - Keramik (99,7%)	hochreine Al <sub>2</sub> O <sub>3</sub> - Keramik (99,7%)	Hastelloy	FEP	PTFE
-----------------------------------	--------------------------	------------------	------------------	---	---	-----------	-----	------

**digital  
burette**

PTFE Compound Halar (ECTFE)	Borosilicate glass 3.3	Halar (ECTFE)	Halar (ECTFE)	highly pure Al <sub>2</sub> O <sub>3</sub> -ceramic (99,7%)	highly pure Al <sub>2</sub> O <sub>3</sub> -ceramic (99,7%)	Hastelloy	FEP	PTFE
-----------------------------------	---------------------------	------------------	------------------	---	---	-----------	-----	------

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for unit ..... (Serial number)

Medium used: .....

Blood, bodily fluids, patholog. samples:

**YES**

**NO**

**Please define in case of Yes**

\_\_\_\_\_

Other hazardous biologic substances:

\_\_\_\_\_

Biologically-degradable materials which can pose a risk:

\_\_\_\_\_

Chemicals containing substances harmful to health:

\_\_\_\_\_

Radioactive materials:

\_\_\_\_\_

Other hazardous and toxic substances:

\_\_\_\_\_

### We hereby declare that

- the unit has been carefully cleaned and decontaminated prior to shipment
- no risks resulting from contamination exist
- we are aware that liability claims can be made against us in the event of damage being caused by contaminated units.

### Address of sender:

Name: \_\_\_\_\_

Date: \_\_\_\_\_

Signature: \_\_\_\_\_